

兔抗 CATSPER3 多克隆抗体

中文名称: 兔抗 CATSPER3 多克隆抗体

英文名称: Anti-CATSPER3 rabbit polyclonal antibody

别 名: CACRC

相关类别: 一抗

储 存: 冷冻(-20℃)

抗 原: CATSPER3

宿 主: Rabbit

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	CatSper3, is a protein which in humans is encoded
	by the CATSPER3 gene. CatSper3 is a member of t
	hecation channels of sperm family of proteins. The
	four proteins in this family together form a Ca2+-p
	ermeant ion channel specific essential for the corre
	ct function of sperm cells. Voltage-gated calcium c
	hannel that plays a central role in calcium-depende
	nt physiological responses essential for successful fe
	rtilization, such as sperm hyperactivation, acrosome
	reaction and chemotaxis towards the oocyte.
Applications:	ELISA, IHC



Name of antibody:	CATSPER3
Immunogen:	Fusion protein of human CATSPER3
Full name:	cation channel, sperm associated 3
Synonyms:	CACRC
SwissProt:	Q86XQ3
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human liver cancer and human thyroid cancer
IHC Recommend dilution:	100-300

