

## 兔抗 DLG4 多克隆抗体

中文名称：兔抗 DLG4 多克隆抗体

英文名称：Anti-DLG4 rabbit polyclonal antibody

别名：PSD95, SAP90, SAP-90

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：DLG4

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the membrane-associated guanylate kinase (MAGUK) family. It heteromultimerizes with another MAGUK protein, DLG2, and is recruited into NMDA receptor and potassium channel clusters. These two MAGUK proteins may interact at postsynaptic sites to form a multimeric scaffold for the clustering of receptors, ion channels, and associated signaling proteins. Multiple transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	DLG4

<b>Immunogen:</b>	Synthetic peptide of human DLG4
<b>Full name:</b>	discs large MAGUK scaffold protein 4
<b>Synonyms:</b>	PSD95; SAP90; SAP-90
<b>SwissProt:</b>	P78352
<b>ELISA Recommended dilution:</b>	2000-10000
<b>IHC positive control:</b>	Human ovarian cancer
<b>IHC Recommend dilution:</b>	50-200

