

兔抗 DLGAP4 多克隆抗体

中文名称：兔抗 DLGAP4 多克隆抗体

英文名称：Anti-DLGAP4 rabbit polyclonal antibody

别名：DLG associated protein 4; DAP4; DLP4; DAP-4; SAPAP4; SAPAP-4

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

抗原：DLGAP4

技术规格

Background:

The product of this gene is a membrane-associated guanylate kinase found at the postsynaptic density in neuronal cells. It is a signaling molecule that can interact with potassium channels and receptors, as well as other signaling molecules. The protein encoded by this gene can interact with PSD-95 through its guanylate kinase domain and may be involved in clustering PSD-95 in the postsynaptic density region. The encoded protein is one of at least four similar proteins that have been found. Alternatively spliced transcript variants encoding different isoforms have been found for this

	gene.
Applications:	ELISA, WB
Name of antibody:	DLGAP4
Immunogen:	Synthetic peptide of human DLGAP4
Full name:	DLG associated protein 4
Synonyms:	DAP4; DLP4; DAP-4; SAPAP4; SAPAP-4
SwissProt:	Q9Y2H0
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	108 kDa
WB Positive control:	A172 cell lysate
WB Recommended dilution:	200-1000

