

兔抗 CAV3 多克隆抗体

中文名称：兔抗 CAV3 多克隆抗体

英文名称：Anti-CAV3 rabbit polyclonal antibody

别名：LQT9; VIP21; LGMD1C; VIP-21

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CAV3

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a caveolin family member, which functions as a component of the caveolae plasma membranes found in most cell types. Caveolin proteins are proposed to be scaffolding proteins for organizing and concentrating certain caveolin-interacting molecules. Mutations identified in this gene lead to interference with protein oligomerization or intra-cellular routing, disrupting caveolae formation and resulting in Limb-Girdle muscular dystrophy type-1C (LGMD-1C), hyperCKemia or rippling muscle disease (RMD). Alternative splicing has been identified for this locus, with
--------------------	--

	inclusion or exclusion of a differentially spliced intron . In addition, transcripts utilize multiple polyA sites and contain two potential translation initiation sites.
Applications:	ELISA, WB, IHC
Name of antibody:	CAV3
Immunogen:	Synthetic peptide of human CAV3
Full name:	caveolin 3
Synonyms:	LQT9; VIP21; LGMD1C; VIP-21
SwissProt:	P56539
IHC positive control:	Human gastric cancer and human thyroid cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	17 kDa
WB Positive control:	Human fetal muscle tissue
WB Recommended dilution:	500-2000



