

6×HIS 抗原

中文名称： 6×HIS 抗原

英文名称： 6×HIS Antigen

相关类别： 抗原

储 存： 冷冻（-20℃）

概述

Synthetic peptide: HHHHHH

技术规格

Full name	HHHHHH Tag
Synonyms	6 His epitope tag
Background	The 6 X histidine motif is used as a tag on many recombinant proteins to facilitate purification.It can then be purified using immobilized metal ions,traditionally nickel,in a procedure termed immobilized-metal affinity chromatography.Further analysis of the protein can be achieved using an antibody specific to the 6 X histidine in techniques such as western blotting.The antibody recognizes the His-tag fused to targeted proteins in transfected or transformed cells.