

## 10X 抗体稀释液(ApplytoWB)

中文名称： 10X 抗体稀释液(Apply to WB)

英文名称： 10X Antibody Diluent for Western Blotting

相关类别： 缓冲液

储 存： 冷藏（2-8℃） 不可冻存

### 概述

Antibody diluent is a 10x concentrated buffer. This buffer is specifically formulated and tested for western blotting applications. Antibody diluent is for diluting both primary and secondary antibodies. The buffer of this antibody diluent is TBS, pH 7.6.

### 使用说明

For 10 ml: 1 ml of antibody diluent (10x) + 9 ml ultra pure water + 5  $\mu$ l Tween20. Now with the current.