

兔抗 CBLN4 多克隆抗体

中文名称：兔抗 CBLN4 多克隆抗体

英文名称：Anti-CBLN4 rabbit polyclonal antibody

别名：CBLNL1

抗原：CBLN4

储存：冷冻（-20℃）避光

宿主：Rabbit

相关类别：一抗

反应种属：Human Mouse

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	This gene encodes a member of a family of small secreted proteins containing C1Q domains. Members of this family are involved in regulation of neuromuscular synapse development. The mouse homolog of the protein encoded by this gene competes with netrin to bind to the deleted in colorectal cancer receptor.
Applications:	WB, IHC
Name of antibody:	CBLN4
Immunogen:	Synthesized peptide derived from internal of human CBLN4.
Full name:	cerebellin 4 precursor

Synonyms :	CBLNL1
SwissProt:	Q9NTU7
IHC positive control:	Human thymus glan tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	22 kDa
WB Positive control:	COS cells lysate
WB Recommended dilution:	500-3000



