



兔抗 CCM2 多克隆抗体

中文名称：兔抗 CCM2 多克隆抗体

英文名称：Anti-CCM2 rabbit polyclonal antibody

别 名：CCM2 scaffolding protein; OSM; C7orf22; PP10187

相关类别：一抗

储 存：冷冻 (-20°C)

宿 主：Rabbit

抗 原：CCM2

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

This gene encodes a scaffold protein that functions in the stress-activated p38 Mitogen-activated protein kinase (MAPK) signaling cascade. The protein interacts with SMAD specific E3 ubiquitin protein ligase 1 (also known as SMURF1) via a phosphotyrosine binding domain to promote RhoA degradation. The protein is required for normal cytoskeletal structure, cell-cell interactions, and lumen formation in endothelial cells. Mutations in this gene result in cerebral cavernous malformations. Multiple transcript variants encoding different isoforms have been found for this gene.

Applications:	ELISA, IHC
Name of antibody:	CCM2
Immunogen:	Fusion protein of human CCM2
Full name:	CCM2 scaffolding protein
Synonyms:	OSM; C7orf22; PP10187
SwissProt:	Q9BSQ5
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	20-100

