

## 兔抗 DNAJC9 多克隆抗体

中文名称：兔抗 DNAJC9 多克隆抗体

英文名称：Anti-DNAJC9 rabbit polyclonal antibody

别名：C3orf41

相关类别：一抗

抗原：DNAJC9

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

|                                    |                                                                                        |
|------------------------------------|----------------------------------------------------------------------------------------|
| <b>Background:</b>                 | May play a role as co-chaperone of the Hsp70 family proteins HSPA1A, HSPA1B and HSPA8. |
| <b>Applications:</b>               | ELISA, WB, IHC                                                                         |
| <b>Name of antibody:</b>           | DNAJC9                                                                                 |
| <b>Immunogen:</b>                  | Synthetic peptide of human DNAJC9                                                      |
| <b>Full name:</b>                  | DnaJ heat shock protein family (Hsp40) member C9                                       |
| <b>Synonyms:</b>                   | JDD1; SB73; HDJC9                                                                      |
| <b>SwissProt:</b>                  | Q8WXX5                                                                                 |
| <b>ELISA Recommended dilution:</b> | 5000-10000                                                                             |
| <b>IHC positive control:</b>       | Human colorectal cancer and Human gastric cancer                                       |
| <b>IHC Recommend dilution:</b>     | 20-100                                                                                 |

|                                 |                            |
|---------------------------------|----------------------------|
| <b>WB Predicted band size:</b>  | 30 kDa                     |
| <b>WB Positive control:</b>     | Human tonsil tissue lysate |
| <b>WB Recommended dilution:</b> | 200-1000                   |



