

兔抗 CCNYL1 多克隆抗体

中文名称：兔抗 CCNYL1 多克隆抗体

英文名称：Anti-CCNYL1 rabbit polyclonal antibody

抗 原：CCNYL1

别 名：cyclin Y like 1;

相关类别：一抗

储 存：冷冻（-20℃）

宿 主：Rabbit

反应种属：Human

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	CCNYL1 belongs to the cyclin family, Cyclin Y subfamily and contains 1 cyclin N terminal domain. The specific function of CCNYL1 is not yet known.
Applications:	ELISA, IHC
Name of antibody:	CCNYL1
Immunogen:	Synthetic peptide of human CCNYL1
Full name:	cyclin Y like 1
SwissProt:	Q8N7R7
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and human lung cancer
IHC Recommend dilution:	25-10



