

兔抗 DTX1 多克隆抗体

中文名称：兔抗 DTX1 多克隆抗体

英文名称： Anti-DTX1 rabbit polyclonal antibody

别 名： hDx-1

抗 原： DTX1

储 存： 冷冻（-20℃）

宿 主： Rabbit

反应种属： Human, Mouse

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Studies in Drosophila have identified this gene as encoding a positive regulator of the Notch-signaling pathway. The human gene encodes a protein of unknown function ; however, it may play a role in basic helix-loop-helix transcription factor activity.
Applications:	ELISA, WB, IHC
Name of antibody:	DTX1
Immunogen:	Synthetic peptide of human DTX1
Full name:	deltex 1, E3 ubiquitin ligase
Synonyms :	hDx-1
SwissProt:	Q86Y01

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer and human gastric cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	67 kDa
WB Positive control:	Human fetal brain tissue and Hela cell
WB Recommended dilution:	200-1000





