

兔抗 CD276 多克隆抗体

中文名称: 兔抗 CD276 多克隆抗体

英文名称: Anti-CD276 rabbit polyclonal antibody

别名: B7H3; B7-H3; B7RP-2; 4Ig-B7-H3

抗原: CD276

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:

The protein encoded by this gene belongs to the immunoglobulin superfamily, and thought to participate in the regulation of T-cell-mediated immune response. Studies show that while the transcript of this gene is ubiquitously expressed in normal tissues and solid tumors, the protein is preferentially expressed only in tumor tissues. Additionally, it was observed that the 3' UTR of this transcript contains a target site for miR29 micro RNA, and there is an inverse correlation between the expression of this protein and miR29 levels, suggesting regulation of expression of this gene product by miR2

	9. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, WB, IHC
Name of antibody:	CD276
Immunogen:	Synthetic peptide of human CD276
Full name:	CD276 molecule
Synonyms:	B7H3; B7-H3; B7RP-2; 4Ig-B7-H3
SwissProt:	Q5ZPR3
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	57 kDa
WB Positive control:	293T , LNCAP , 231 and PC-3 cell lysates
WB Recommended dilution:	200-1000



