

兔抗 DUS2 多克隆抗体

中文名称：兔抗 DUS2 多克隆抗体

英文名称：Anti-DUS2 rabbit polyclonal antibody

别名：dihydrouridine synthase 2; SMM1; DUS2L; URLC8

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：DUS2

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a cytoplasmic protein that catalyzes the conversion of uridine residues to dihydrouridine in the D-loop of tRNA. The resulting modified bases confer enhanced regional flexibility to tRNA. The encoded protein may increase the rate of translation by inhibiting an interferon-induced protein kinase. This gene has been implicated in pulmonary carcinogenesis. Alternatively spliced transcript variants have been described for this gene.
Applications:	ELISA, WB, IHC
Name of antibody:	DUS2

Immunogen:	Fusion protein of human DUS2
Full name:	dihydrouridine synthase 2
Synonyms:	SMM1; DUS2L; URLC8
SwissProt:	Q9NX74
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	30-150
WB Predicted band size:	55 kDa
WB Positive control:	HepG2 cell and Human Bladder transitional cell carcinoma grade 2-3 tissue lysates
WB Recommended dilution:	200-1000



