

兔抗 DYNC1LI2 多克隆抗体

中文名称：兔抗 DYNC1LI2 多克隆抗体

英文名称：Anti-DYNC1LI2 rabbit polyclonal antibody

别名：LIC2; DNCL12

相关类别：一抗

储存：冷冻（-20℃）

抗原：DYNC1LI2

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Acts as one of several non-catalytic accessory components of the cytoplasmic dynein 1 complex that are thought to be involved in linking dynein to cargos and to adapter proteins that regulate dynein function. Cytoplasmic dynein 1 acts as a motor for the intracellular retrograde motility of vesicles and organelles along microtubules. May play a role in binding dynein to membranous organelles or chromosomes.
Applications:	ELISA, WB, IHC
Name of antibody:	DYNC1LI2
Immunogen:	Synthetic peptide of human DYNC1LI2

Full name:	dynein cytoplasmic 1 light intermediate chain 2
Synonyms:	LIC2; DNCLI2
SwissProt:	O43237
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer and Human prostate cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	54 kDa
WB Positive control:	Human fetal brain tissue and Mouse brain tissue, PC-3 cell lysates
WB Recommended dilution:	200-1000



