

## 兔抗 DZIP1 多克隆抗体

中文名称: 兔抗 DZIP1 多克隆抗体

英文名称: Anti-DZIP1 rabbit polyclonal antibody

别 名: DZIP; DZIPt1

相关类别: 一抗

储 存: 冷冻 (-20℃)

抗 原: DZIP1

宿 主: Rabbit

反应种属: Human

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

## 技术规格

**Background:** 

DZIP1(DAZ interacting protein 1) is also known as DZIP or DZIP2 and is a 867 amino acid protein which is expressed as 3 isoforms, designated DZIPb, DZIPt1 and DZIPt2. DZIP 1 is localized to testis, oocytes, ovary and fetal brain, as well as in embryonic stem cells and germ cells. In testis, DZIP1 is localized to the nucleus and also shows some cy toplasmic distribution in spermatogonia. DZIP1 belongs to the C2H2-type zinc-finger protein family, and, characterists of the C2H2-type zinc-finger protein family, DZIP1 contain s one C2H2-type zinc-finger region through which it is thought to interact with DAZ, an interaction that promotes



	spermiogenesis. DZIP1 expression is not found in those af flicted with Sertoli cell-only syndrome (characterized by th e absence of germ cells in the testis), suggesting that the lack of DZIP1 may be involved in the pathogenesis of Ser toli cell-only syndrome-induced male sterility
Applications:	ELISA, IHC
Name of antibody:	DZIP1
Immunogen:	Synthetic peptide of human DZIP1
Full name:	DAZ interacting zinc finger protein 1
Synonyms:	DZIP; DZIPt1
SwissProt:	Q86YF9
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer;
IHC Recommend dilution:	25-100

