

兔抗 CEACAM8 多克隆抗体

中文名称：兔抗 CEACAM8 多克隆抗体

英文名称：Anti-CEACAM8 rabbit polyclonal antibody

别名：CD67; CGM6; CD66b; NCA-95

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CEACAM8

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Cell surface glycoprotein that plays a role in cell adhesion in a calcium-independent manner (PubMed:8776764, PubMed:2022629, PubMed:11590190). Mediates heterophilic cell adhesion with other carcinoembryonic antigen-related cell adhesion molecules, such as CEACAM6 (PubMed:8776764, PubMed:2022629, PubMed:11590190). Heterophilic interaction with CEACAM8 occurs in activated neutrophils (PubMed:8776764).
Applications:	ELISA, IHC
Name of antibody:	CEACAM8
Immunogen:	Fusion protein of human CEACAM8

Full name:	CEA cell adhesion molecule 8
Synonyms:	CD67; CGM6; CD66b; NCA-95
SwissProt:	P31997
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and Human tonsil
IHC Recommend dilution:	500-1000



