

兔抗 EGFR(Phospho-Ser695) 多克隆抗体

中文名称：兔抗 EGFR(Phospho-Ser695) 多克隆抗体

英文名称：Anti-EGFR(Phospho-Ser695) rabbit polyclonal antibody

别名：ERBB; HER1; mENA; ERBB1; PIG61

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：EGFR(Phospho-Ser695)

反应种属：Human Mouse Rat

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	The protein encoded by this gene is a transmembrane glycoprotein that is a member of the protein kinase superfamily. This protein is a receptor for members of the epidermal growth factor family. EGFR is a cell surface protein that binds to epidermal growth factor. Binding of the protein to a ligand induces receptor dimerization and tyrosine autophosphorylation and leads to cell proliferation. Mutations in this gene are associated with lung cancer.
Applications:	WB
Name of antibody:	EGFR(Phospho-Ser695)
Immunogen:	Peptide sequence around phosphorylation site of Serine 695

	(T-P-S(p)-G-E) derived from Human EGFR.
Full name:	epidermal growth factor receptor
Synonyms :	ERBB; HER1; mENA; ERBB1; PIG61
SwissProt:	P00533
WB Predicted band size:	134 kDa
WB Positive control:	293 cells lysates treated with EGF
WB Recommended dilution:	500-1000

