

## 兔抗 EGFR(Phospho-Ser1071) 多克隆抗体

中文名称：兔抗 EGFR(Phospho-Ser1071) 多克隆抗体

英文名称：Anti-EGFR(Phospho-Ser1071) rabbit polyclonal antibody

别名：ERBB; HER1; mENA; ERBB1; PIG61

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：EGFR(Phospho-Ser1071)

反应种属：Human Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a transmembrane glycoprotein that is a member of the protein kinase superfamily. This protein is a receptor for members of the epidermal growth factor family. EGFR is a cell surface protein that binds to epidermal growth factor. Binding of the protein to a ligand induces receptor dimerization and tyrosine autophosphorylation and leads to cell proliferation. Mutations in this gene are associated with lung cancer.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	EGFR(Phospho-Ser1071)
<b>Immunogen:</b>	Peptide sequence around phosphorylation site of Serine 10

	71 (Y-S-S(p)-D-P) derived from Human EGFR.
<b>Full name:</b>	epidermal growth factor receptor
<b>Synonyms :</b>	ERBB; HER1; mENA; ERBB1; PIG61
<b>SwissProt:</b>	P00533
<b>WB Predicted band size:</b>	134 kDa
<b>WB Positive control:</b>	Rat renin lysates
<b>WB Recommended dilution:</b>	500-1000

