

兔抗 CELA3A 多克隆抗体

中文名称：兔抗 CELA3A 多克隆抗体

英文名称：Anti-CELA3A rabbit polyclonal antibody

别 名：chymotrypsin like elastase family member 3A; ELA3; ELA3A

相关类别：一抗

储 存：冷冻（-20℃）

宿 主：Rabbit

抗 原：CELA3A

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Elastases form a subfamily of serine proteases that hydrolyze many proteins in addition to elastin. Humans have six elastase genes which encode the structurally similar proteins elastase 1, 2, 2A, 2B, 3A, and 3B. Unlike other elastases, elastase 3A has little elastolytic activity. Like most of the human elastases, elastase 3A is secreted from the pancreas as a zymogen and, like other serine proteases such as trypsin, chymotrypsin and kallikrein, it has a digestive function in the intestine. Elastase 3A preferentially cleaves proteins after alanine residues. Elastase 3A may also function in the intestinal transport
--------------------	---

	t and metabolism of cholesterol. Both elastase 3A and elastase 3B have been referred to as protease E and as elastase 1.
Applications:	ELISA, IHC
Name of antibody:	CELA3A
Immunogen:	Fusion protein of human CELA3A
Full name:	chymotrypsin like elastase family member 3A
Synonyms:	ELA3; ELA3A
SwissProt:	P09093
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	20-100

