

兔抗 EHMT2 多克隆抗体

中文名称：兔抗 EHMT2 多克隆抗体

英文名称：Anti-EHMT2 rabbit polyclonal antibody

别名：euchromatic histone lysine methyltransferase 2; G9A; BAT8; GAT8; NG36; KMT1C; C6orf30

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：EHMT2

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a methyltransferase that methylates lysine residues of histone H3. Methylation of H3 at lysine 9 by this protein results in recruitment of additional epigenetic regulators and repression of transcription. This gene was initially thought to be two different genes, NG36 and G9a, adjacent to each other in the HLA locus. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	EHMT2
Immunogen:	Synthetic peptide of human EHMT2

Full name:	euchromatic histone lysine methyltransferase 2
Synonyms:	G9A; BAT8; GAT8; NG36; KMT1C; C6orf30
SwissProt:	Q96KQ7
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human ovarian cancer and Human cervical cancer
IHC Recommend dilution:	25-100



