

兔抗 EIF2B2 多克隆抗体

中文名称：兔抗 EIF2B2 多克隆抗体

英文名称：Anti-EIF2B2 rabbit polyclonal antibody

别名：eukaryotic translation initiation factor 2B subunit beta; EIF2B; EIF-2Bbeta

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：EIF2B2

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes the beta subunit of eukaryotic initiation factor-2B (EIF2B). EIF2B is involved in protein synthesis and exchanges GDP and GTP for its activation and deactivation.
Applications:	ELISA, IHC
Name of antibody:	EIF2B2
Immunogen:	Fusion protein of human EIF2B2
Full name:	eukaryotic translation initiation factor 2B subunit beta
Synonyms:	EIF2B; EIF-2Bbeta
SwissProt:	P49770

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human gastric cancer and Human prostate cancer
IHC Recommend dilution:	50-200



