

## 兔抗 CFAP45 多克隆抗体

中文名称：兔抗 CFAP45 多克隆抗体

英文名称：Anti-CFAP45 rabbit polyclonal antibody

别名：NESG1; CCDC19

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CFAP45

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

#### Background:

CCDC19 is a 466 amino acid protein encoded by a gene mapping to human chromosome 1. Chromosome 1 is the largest human chromosome, spanning about 260 million base pairs and making up 8% of the human genome. There are about 3,000 genes on chromosome 1 and, considering the great number of genes, there are also a large number of diseases associated with chromosome 1. Notably, the rare aging disease Hutchinson-Gilford progeria is associated with the LMNA gene which encodes Lamin A. When defective, the LMNA gene product can build up in the nucleus and cause characteristic nuclear abnormalities.

	lear blebs. The mechanism of rapidly enhanced aging is unclear and is a topic of continuing exploration. The M UTYH gene is located on chromosome 1 and is partially responsible for familial adenomatous polyposis. Stickler syndrome, Parkinson's, Gaucher disease and Usher syndrome are also associated with chromosome 1.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CFAP45
<b>Immunogen:</b>	Synthetic peptide of humanCFAP45
<b>Full name:</b>	cilia and flagella associated protein 45
<b>Synonyms :</b>	NESG1; CCDC19
<b>SwissProt:</b>	Q9UL16
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	25-100

