

兔抗 CFAP298 多克隆抗体

中文名称：兔抗 CFAP298 多克隆抗体

英文名称：Anti-CFAP298 rabbit polyclonal antibody

别名：cilia and flagella associated protein 298; Kur; FB

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CFAP298

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a protein that plays a critical role in dynein arm assembly and motile cilia function. Mutations in this gene result in primary ciliary dyskinesia. Naturally occurring readthrough transcription occurs from this locus to the downstream t-complex 10 like (TCP10L) gene. [provided by RefSeq, Apr 2017]
Applications:	ELISA, WB, IHC
Name of antibody:	CFAP298
Immunogen:	Fusion protein of human CFAP298
Full name:	cilia and flagella associated protein 298
Synonyms:	Kur; FBB18; CILD26; C21orf48; C21orf59

SwissProt:	P57076
IHC positive control:	Human prostate cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	33 kDa
WB Positive control:	Human testis tissue lysate
WB Recommended dilution:	200-1000



