

兔抗 EIF5 多克隆抗体

中文名称：兔抗 EIF5 多克隆抗体

英文名称：Anti-EIF5 rabbit polyclonal antibody

别名：EIF-5; EIF-5A

相关类别：一抗

抗原：EIF5

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Eukaryotic translation initiation factor-5 (EIF5) interacts with the 40S initiation complex to promote hydrolysis of bound GTP with concomitant joining of the 60S ribosomal subunit to the 40S initiation complex. The resulting functional 80S ribosomal initiation complex is then active in peptidyl transfer and chain elongations
Applications:	ELISA, IHC
Name of antibody:	EIF5
Immunogen:	Fusion protein of human EIF5
Full name:	eukaryotic translation initiation factor 5

Synonyms :	EIF-5; EIF-5A
SwissProt:	P55010
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human breast cancer and human ovarian cancer
IHC Recommend dilution:	25-100

