

## 兔抗 ELANE 多克隆抗体

- 中文名称: 兔抗 ELANE 多克隆抗体
- 英文名称: Anti-ELANE rabbit polyclonal antibody
- 别 名: elastase, neutrophil expressed; GE; NE; HLE; HNE; ELA2; SCN1; PMN-E
- 相关类别: 一抗
- 抗 原: ELANE
- 储 存: 冷冻 (-20℃)
- 宿 主: Rabbit
- 反应种属: Human
- 标记物: Unconjugate
- 克隆类型: rabbit polyclonal

## 技术规格

Background:	Elastases form a subfamily of serine proteases that hydrol yze many proteins in addition to elastin. Humans have si x elastase genes which encode structurally similar protein s. The encoded preproprotein is proteolytically processed to generate the active protease. Following activation, this protease hydrolyzes proteins within specialized neutrophil lysosomes, called azurophil granules, as well as proteins
	lysosomes, called azurophil granules, as well as proteins of the extracellular matrix. The enzyme may play a role i n degenerative and inflammatory diseases through proteo lysis of collagen-IV and elastin. This protein also degrade s the outer membrane protein A (OmpA) of E. coli as we



	Il as the virulence factors of such bacteria as Shigella, Sal monella and Yersinia. Mutations in this gene are associat ed with cyclic neutropenia and severe congenital neutrop enia (SCN). This gene is present in a gene cluster on chr omosome 19.
Applications:	ELISA, IHC
Name of antibody:	ELANE
Immunogen:	Synthetic peptide of human ELANE
Full name:	elastase, neutrophil expressed
Synonyms:	GE; NE; HLE; HNE; ELA2; SCN1; PMN-E
SwissProt:	P08246
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer and Human prostate cancer
IHC Recommend dilution:	30-150





