

## 兔抗 STMN2/STMN3/STMN4 多克隆抗体

中文名称：兔抗 STMN2/STMN3/STMN4 多克隆抗体

英文名称：Anti-STMN2/STMN3/STMN4 rabbit polyclonal antibody

别名：SCG10, SCGN10; SCLIP; RB3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：STMN2/STMN3/STMN4

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the stathmin family of phosphoproteins. Stathmin proteins function in microtubule dynamics and signal transduction. The encoded protein plays a regulatory role in neuronal growth and is also thought to be involved in osteogenesis. Reductions in the expression of this gene have been associated with Down's syndrome and Alzheimer's disease. Alternatively spliced transcript variants have been observed for this gene. A pseudogene of this gene is located on the long arm of chromosome 6.
<b>Applications:</b>	ELISA, WB, IHC

<b>Name of antibody:</b>	STMN2/STMN3/STMN4
<b>Immunogen:</b>	Synthetic peptide of human STMN2/STMN3/STMN4
<b>Full name:</b>	stathmin-like 2/2/3
<b>Synonyms :</b>	SCG10, SCGN10; SCLIP; RB3
<b>SwissProt:</b>	Q93045/Q9NZ72/Q9H169
<b>ELISA Recommended dilution:</b>	1000-5000
<b>IHC positive control:</b>	Human thyroid cancer and Human lung cancer
<b>IHC Recommend dilution:</b>	50-200
<b>WB Predicted band size:</b>	21 kDa
<b>WB Positive control:</b>	Human liver cancer tissue
<b>WB Recommended dilution:</b>	500-2000



