

兔抗 CDC20 多克隆抗体

中文名称：兔抗 CDC20 多克隆抗体

英文名称： Anti-CDC20 rabbit polyclonal antibody

别名： cell division cycle 20; CDC20A; p55CDC; bA276H19.3

抗原： CDC20

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	CDC20 appears to act as a regulatory protein interacting with several other proteins at multiple points in the cell cycle. It is required for two microtubule-dependent processes, nuclear movement prior to anaphase and chromosome separation.
Applications:	ELISA, IHC
Name of antibody:	CDC20
Immunogen:	Fusion protein of human CDC20
Full name:	cell division cycle 20
Synonyms:	CDC20A; p55CDC; bA276H19.3
SwissProt:	Q12834

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	70-350

