

兔抗 EN1 多克隆抗体

中文名称：兔抗 EN1 多克隆抗体

英文名称：Anti-EN1 rabbit polyclonal antibody

别名：engrailed homeobox 1

抗原：EN1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Homeobox-containing genes are thought to have a role in controlling development. In Drosophila, the 'engrailed' (en) gene plays an important role during development in segmentation, where it is required for the formation of posterior compartments. Different mutations in the mouse homologs, En1 and En2, produced different developmental defects that frequently are lethal. The human engrailed homologs 1 and 2 encode homeodomain-containing proteins and have been implicated in the control of pattern formation during development of the central nervous system.
--------------------	---

Applications:	ELISA, IHC
Name of antibody:	EN1
Immunogen:	Synthetic peptide of human EN1
Full name:	engrailed homeobox 1
SwissProt:	Q05925
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	20-100

