

兔抗 ENDOG 多克隆抗体

中文名称：兔抗 ENDOG 多克隆抗体

英文名称：Anti-ENDOG rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻（-20℃）

宿 主：Rabbit

抗 原：ENDOG

反应种属：Human

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a nuclear encoded endonuclease that is localized in the mitochondrion. The encoded protein is widely distributed among animals and cleaves DNA at GC tracts. This protein is capable of generating the RNA primers required by DNA polymerase gamma to initiate replication of mitochondrial DNA.
Applications:	ELISA, IHC
Name of antibody:	ENDOG
Immunogen:	Synthetic peptide of human ENDOG
Full name:	endonuclease G
SwissProt:	Q14249
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human thyroid cancer and Human brain

IHC Recommend dilution:	25-100
-------------------------	--------

