

兔抗 ENPEP 多克隆抗体

中文名称：兔抗 ENPEP 多克隆抗体

英文名称：Anti-ENPEP rabbit polyclonal antibody

别名：APA; CD249; gp160

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ENPEP

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Appears to have a role in the catabolic pathway of the renin-angiotensin system. Probably plays a role in regulating growth and differentiation of early B-lineage cells.
Applications:	ELISA, IHC
Name of antibody:	ENPEP
Immunogen:	Synthetic peptide of human ENPEP
Full name:	glutamyl aminopeptidase
Synonyms:	APA; CD249; gp160
SwissProt:	Q07075
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer

IHC Recommend dilution:	50-300
-------------------------	--------

