

PAK4 抗原(重组蛋白)

中文名称: PAK4 抗原(重组蛋白)

英文名称: PAK4 Antigen (Recombinant Protein)

储 存: 冷冻(-20℃)

相关类别: 抗原

概述

Fusion protein corresponding to a region derived from 321-572 amino acids of human PAK4

技术规格

Full name:	p21 protein (Cdc42/Rac)-activated kinase 4
Swissprot:	O96013
Gene Accession:	BC011368
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	PAK proteins, a family of serine/threonine p21-activating kinases, include PAK1, PAK2, PAK3 and PAK4. PAK proteins are critical ef fectors that link Rho GTPases to cytoskeleton reorganization and nuclear signaling. They serve as targets for the small GTP bindin g proteins Cdc42 and Rac and have been implicated in a wide r ange of biological activities. PAK4 interacts specifically with the GTP-bound form of Cdc42Hs and weakly activates the JNK famil y of MAP kinases. PAK4 is a mediator of filopodia formation and may play a role in the reorganization of the actin cytoskeleto n. Multiple alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.

