

兔抗 EPGN 多克隆抗体

中文名称：兔抗 EPGN 多克隆抗体

英文名称：Anti-EPGN rabbit polyclonal antibody

别名：EPG; PRO9904; ALGV3072

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：EPGN

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the epidermal growth factor family. Members of this family are ligands for the epidermal growth factor receptor and play a role in cell survival, proliferation and migration. This protein has been reported to have high mitogenic activity but low affinity for its receptor. Expression of this transcript and protein have been reported in cancer specimens of the breast, bladder, and prostate. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, WB
Name of antibody:	EPGN

Immunogen:	Synthetic peptide of human EPGN
Full name:	epithelial mitogen
Synonyms :	EPG; PRO9904; ALGV3072
SwissProt:	Q6UW88
ELISA Recommended dilution:	1000-2000
WB Predicted band size:	17 kDa
WB Positive control:	Jurkat, 293T and hela cells
WB Recommended dilution:	200-1000



