

## POU2F1 抗原（重组蛋白）

中文名称：POU2F1 抗原（重组蛋白）

英文名称：POU2F1 Antigen (Recombinant Protein)

别 名：OCT1, OTF1, oct-1B

储 存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 1-290 amino acids of human POU2F1

技术规格

<b>full name:</b>	POU class 2 homeobox 1
<b>Synonyms:</b>	OCT1, OTF1, oct-1B
<b>Swissprot:</b>	P14859
<b>Gene Accession:</b>	BC001664
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The OCT1 transcription factor was among the first identified members of the POU transcription factor family. Members of this family contain the POU domain, a 160-amino acid region necessary for DNA binding to the octameric sequence ATGCAAAT.