

PRMT7 抗原(重组蛋白)

中文名称: PRMT7 抗原(重组蛋白)

英文名称: PRMT7 Antigen (Recombinant Protein)

储 存: 冷冻(-20℃)

相关类别: 抗原

概述

Fusion protein corresponding to a region derived from 526-692 amino acids of human PRMT7

技术规格

Full name:	Protein arginine methyltransferase 7
Swissprot:	Q9NVM4
Gene Accession:	BC000146
Purity:	>85%, as determined by Coomassie blue stained SDS-PA GE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Protein arginine methyltransferase 7 is a protein that in humans is encoded by the PRMT7 gene. Arginine methyl ation is an apparently irreversible protein modification ca talyzed by arginine methyltransferases, such as PMT7, usi ng S-adenosylmethionine (AdoMet) as the methyl donor. Arginine methylation is implicated in signal transduction, RNA transport, and RNA splicing. May be involved in eto poside-induced cytotoxicity, a chemotherapeutic agent fre quently used for testicular cancer and small-cell lung can cer that can cause cytotoxicity in the treatment of other



cancers. Down-regulation confers increased sensitivity to
the Top1 inhibitor camptothecin (CPT).