

PSMD3 抗原(重组蛋白)

中文名称: PSMD3 抗原(重组蛋白)

英文名称: PSMD3 Antigen (Recombinant Protein)

别 名: S3; P58; RPN3; TSTA2

储 存: 冷冻(-20℃)

相关类别: 抗原

概述

Fusion protein corresponding to C terminal 250 amino acids of human PSMD3

技术规格

| Full name: | proteasome (prosome, macropain) 26S subunit, non-ATPase, 3 |
|--------------------|--|
| Synonyms: | S3; P58; RPN3; TSTA2 |
| Swissprot: | O43242 |
| Gene Accession: | BC012302 |
| Purity: | >85%, as determined by Coomassie blue stained SDS-PAGE |
| Expression system: | Escherichia coli |
| Tags: | His tag C-Terminus, GST tag N-Terminus |
| Background: | The 26S proteasome is a multicatalytic proteinase complex with a highly ordered structure composed of 2 complexes, a 20S core a nd a 19S regulator. The 20S core is composed of 4 rings of 28 non-identical subunits; 2 rings are composed of 7 alpha subunits and 2 rings are composed of 7 beta subunits. The 19S regulator is composed of a base, which contains 6 ATPase subunits and 2 non-ATPase subunits, and a lid, which contains up to 10 non-AT |



Pase subunits. Proteasomes are distributed throughout eukaryotic cells at a high concentration and cleave peptides in an ATP/ubiq uitin-dependent process in a non-lysosomal pathway. This gene e ncodes a member of the proteasome subunit S3 family that func tions as one of the non-ATPase subunits of the 19S regulator lid . Single nucleotide polymorphisms in this gene are associated with neutrophil count.