

兔抗 EPHX3 多克隆抗体

中文名称：兔抗 EPHX3 多克隆抗体

英文名称：Anti-EPHX3 rabbit polyclonal antibody

别 名：EH3; ABHD9

抗 原：EPHX3

储 存：冷冻（-20℃）避光

宿 主：Rabbit

相关类别：一抗

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	Belongs to the AB hydrolase superfamily.
Applications:	WB
Name of antibody:	EPHX3
Immunogen:	Synthesized peptide derived from internal of human ABHD9.
Full name:	epoxide hydrolase 3
Synonyms :	EH3; ABHD9
SwissProt:	Q9H6B9
WB Predicted band size:	41 kDa
WB Positive control:	MCF-7 cells lysate
WB Recommended dilution:	500-3000

