

兔抗 EPN2 多克隆抗体

中文名称: 兔抗 EPN2 多克隆抗体

英文名称: Anti-EPN2 rabbit polyclonal antibody

别 名: epsin 2; EHB21

相关类别: 一抗

抗 原: EPN2

储 存: 冷冻(-20℃)

宿 主: Rabbit

反应种属: Human, Rat

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a protein which interacts with clathrin and adaptor-related protein complex 2, alpha 1 subunit. The protein is found in a brain-derived clathrin-coated vesicle fraction and localizes to the peri-Golgi region and the cell periphery. The protein is thought to be involved in clathrin-mediated endocytosis. Alternate splicing of this gene results in multiple transcript variants encoding different isoforms.
Applications:	ELISA, IHC
Name of antibody:	EPN2
Immunogen:	Synthetic peptide of human EPN2



Full name:	epsin 2
Synonyms:	EHB21
SwissProt:	O95208
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	20-100

