

兔抗 EPO 多克隆抗体

中文名称：兔抗 EPO 多克隆抗体

英文名称：Anti-EPO rabbit polyclonal antibody

别名：EP; MVCD2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：EPO

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene is a member of the EPO/TPO family and encodes a secreted, glycosylated cytokine composed of four alpha helical bundles. The protein is found in the plasma and regulates red cell production by promoting erythroid differentiation and initiating hemoglobin synthesis. This protein also has neuroprotective activity against a variety of potential brain injuries and antiapoptotic functions in several tissue types.
Applications:	ELISA, IHC
Name of antibody:	EPO
Immunogen:	Synthetic peptide of human EPO

Full name:	erythropoietin
Synonyms :	EP; MVCD2
SwissProt:	P01588
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human liver cancer and Human thyroid cancer
IHC Recommend dilution:	25-100

