

兔抗 EPS8L2 多克隆抗体

中文名称：兔抗 EPS8L2 多克隆抗体

英文名称： Anti-EPS8L2 rabbit polyclonal antibody

别 名： EPS8R2

储 存： 冷冻（-20℃） 避光

抗 原： EPS8L2

宿 主： Rabbit

反应种属： Human Mouse

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes a member of the EPS8 gene family. The encoded protein, like other members of the family, is thought to link growth factor stimulation to actin organization, generating functional redundancy in the pathways that regulate actin cytoskeletal remodeling.
Applications:	WB
Name of antibody:	EPS8L2
Immunogen:	Synthesized peptide derived from internal of human ES8L2.
Full name:	EPS8-like 2
Synonyms :	EPS8R2
SwissProt:	Q9H6S3

WB Predicted band size:	81 kDa
WB Positive control:	Hela cells, Jurkat cells and Huvec cells lysates
WB Recommended dilution:	500-3000

