

SNAI1 抗原（重组蛋白）

中文名称： SNAI1 抗原（重组蛋白）

英文名称： SNAI1 Antigen (Recombinant Protein)

别 名： SNA; SNAH; SNAIL; SLUGH2; SNAIL1; dJ710H13.1

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 1-264 amino acids of human SNAI1

技术规格

Full name:	snail family transcriptional repressor 1
Synonyms:	SNA; SNAH; SNAIL; SLUGH2; SNAIL1; dJ710H13.1
Swissprot:	O95863
Gene Accession:	BC012910
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The Drosophila embryonic protein snail is a zinc finger transcriptional repressor which downregulates the expression of ectodermal genes within the mesoderm. The nuclear protein encoded by this gene is structurally similar to the Drosophila snail protein, and is also thought to be critical for mesoderm formation in the developing embryo. At least two variants of a similar processed pseudogene have been found on chromosome 2.