

兔抗 ERF 多克隆抗体

中文名称：兔抗 ERF 多克隆抗体

英文名称：Anti-ERF rabbit polyclonal antibody

别名：PE2; CRS4; PE-2

相关类别：一抗

储存：冷冻（-20℃）避光

抗原：ERF

宿主：Rabbit

反应种属：Human Mouse

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	Members of the ETS family of transcription factors, such as ERF, regulate cell proliferation and differentiation. They share a highly conserved DNA-binding domain, the ETS domain, that recognizes the sequence GGAA/T (de Castro et al., 1997 [PubMed 9192842]). For further information on ETS transcription factors, see ETS1 (MIM 164720).
Applications:	IHC
Name of antibody:	ERF
Immunogen:	Synthesized peptide derived from internal of human ERF.
Full name:	Ets2 repressor factor
Synonyms :	PE2; CRS4; PE-2

SwissProt:	P50548
IHC positive control:	Human breast carcinoma tissue
IHC Recommend dilution:	50-100

