

兔抗 ERMP1 多克隆抗体

中文名称：兔抗 ERMP1 多克隆抗体

英文名称：Anti-ERMP1 rabbit polyclonal antibody

别名：FXNA; KIAA1815; bA207C16.3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ERMP1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Within the ovary, required for the organization of somatic cells and oocytes into discrete follicular structures.
Applications:	ELISA, IHC
Name of antibody:	ERMP1
Immunogen:	Fusion protein of human ERMP1
Full name:	endoplasmic reticulum metalloproteinase 1
Synonyms :	FXNA; KIAA1815; bA207C16.3
SwissProt:	Q7Z2K6
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human breast cancer and human thyroid cancer
IHC Recommend dilution:	25-100



