

兔抗 ERP27 多克隆抗体

中文名称：兔抗 ERP27 多克隆抗体

英文名称：Anti-ERP27 rabbit polyclonal antibody

别名：endoplasmic reticulum protein 27; PDIA8; C12orf46

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ERP27

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	ERP27 is a noncatalytic member of the protein disulfide isomerase (PDI; see MIM 608012) family of endoplasmic reticulum (ER) proteins (Alanen et al., 2006 [PubMed 16940051]).
Applications:	ELISA, IHC
Name of antibody:	ERP27
Immunogen:	Fusion protein of human ERP27
Full name:	endoplasmic reticulum protein 27
Synonyms:	PDIA8; C12orf46
SwissProt:	Q96DN0
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human cervical cancer
IHC Recommend dilution:	40-200

