

兔抗 ESM1 多克隆抗体

中文名称：兔抗 ESM1 多克隆抗体

英文名称：Anti-ESM1 rabbit polyclonal antibody

别名：endothelial cell specific molecule 1; endocan

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ESM1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a secreted protein which is mainly expressed in the endothelial cells in human lung and kidney tissues. The expression of this gene is regulated by cytokines, suggesting that it may play a role in endothelium-dependent pathological disorders. The transcript contains multiple polyadenylation and mRNA instability signals. Two transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	ESM1
Immunogen:	Fusion protein of human ESM1

Full name:	endothelial cell specific molecule 1
Synonyms:	endocan
SwissProt:	Q9NQ30
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer and Human colorectal cancer
IHC Recommend dilution:	50-300



