

兔抗 ETS2 多克隆抗体

中文名称：兔抗 ETS2 多克隆抗体

英文名称：Anti-ETS2 rabbit polyclonal antibody

别名：ETS2IT1

相关类别：一抗

储存：冷冻（-20℃）

抗原：ETS2

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a transcription factor which regulates genes involved in development and apoptosis. The encoded protein is also a protooncogene and shown to be involved in regulation of telomerase. A pseudogene of this gene is located on the X chromosome. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, WB
Name of antibody:	ETS2
Immunogen:	Fusion protein of human ETS2
Full name:	v-ets avian erythroblastosis virus E26 oncogene homolog 2

Synonyms :	ETS2IT1
SwissProt:	P15036
ELISA Recommended dilution:	2000-5000
WB Predicted band size:	53 kDa
WB Positive control:	MCF-7 cells
WB Recommended dilution:	500-2000

