

CXCR4 抗原(重组蛋白)

中文名称: CXCR4 抗原(重组蛋白)

英文名称: CXCR4 Antigen (Recombinant Protein)

别 名: FB22; HM89; LAP3; LCR1; NPYR; WHIM; CD184; LESTR; NPY3R; NPYRL; HSY3RR; NPYY3R; D2S201E

- 储存: 冷冻(-20℃)
- 相关类别: 抗原
- 概 述:

Fusion protein corresponding to a region derived from 262-352 amino acids of human CXCR4

技术规格:

Full name:	chemokine (C-X-C motif) receptor 4
Synonyms:	FB22; HM89; LAP3; LCR1; NPYR; WHIM; CD184; LESTR; NPY3R; NPYRL; HSY3RR; NPYY3R; D2S201E
Swissprot:	P61073
Gene Accession:	BC020968
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a CXC chemokine receptor specific for stromal cell -derived factor-1. The protein has 7 transmembrane regions and is lo cated on the cell surface. It acts with the CD4 protein to support HIV entry into cells and is also highly expressed in breast cancer cells. M utations in this gene have been associated with WHIM (warts, hypoga



mmaglobulinemia, infections, and myelokathexis) syndrome. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.