

EDN2 抗原（重组蛋白）

中文名称： EDN2 抗原（重组蛋白）

英文名称： EDN2 Antigen (Recombinant Protein)

别 名： ET2; PPET2

储 存： 冷冻（-20℃）

相关类别： 抗原

概 述：

Fusion protein corresponding to a region derived from 70-178 amino acids of human EDN2

技术规格：

Full name:	endothelin 2
Synonyms:	ET2; PPET2
Swissprot:	P20800
Gene Accession:	BC034393
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the endothelin protein family of secretory vasoconstrictive peptides. The preproprotein is processed to a short mature form which functions as a ligand for the endothelin receptors that initiate intracellular signaling events. This gene product is involved in a wide range of biological processes, such as hypertension and ovulation.