

## FES 抗原(重组蛋白)

- 中文名称: FES 抗原(重组蛋白)
- 英文名称: FES Antigen (Recombinant Protein)
- 别名: FPS
- 储存: 冷冻(-20℃)
- 相关类别: 抗原
- 概 述:

Fusion protein corresponding to C terminal 250 amino acids of human FES

## 技术规格:

Full name:	feline sarcoma oncogene
Synonyms:	FPS
Swissprot:	P07332
Gene Accession:	BC035357
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes the human cellular counterpart of a felin e sarcoma retrovirus protein with transforming capabilities. T he gene product has tyrosine-specific protein kinase activity and that activity is required for maintenance of cellular trans formation. Its chromosomal location has linked it to a specif ic translocation event identified in patients with acute promy elocytic leukemia but it is also involved in normal hematopo



iesis as well as growth factor and cytokine receptor signalin g. Alternative splicing results in multiple variants encoding d ifferent isoforms.