

FKBP8 抗原(重组蛋白)

- 中文名称: FKBP8 抗原(重组蛋白)
- 英文名称: FKBP8 Antigen (Recombinant Protein)
- 别名: FKBP38; FKBPr38
- 储存: 冷冻(-20℃)
- 相关类别: 抗原
- 概 述:

Fusion protein corresponding to a region derived from 76-325 amino acids of human FKBP8

技术规格:

Full name:	FK506 binding protein 8, 38kDa
Synonyms:	FKBP38; FKBPr38
Swissprot:	Q14318
Gene Accession:	BC009966
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is a member of the imm unophilin protein family, which play a role in immunoregulat ion and basic cellular processes involving protein folding an d trafficking. Unlike the other members of the family, this e ncoded protein does not seem to have PPIase/rotamase acti vity. It may have a role in neurons associated with memory function.