

GLIPR1 抗原(重组蛋白)

- 中文名称: GLIPR1 抗原(重组蛋白)
- 英文名称: GLIPR1 Antigen (Recombinant Protein)
- 别名: GLIPR; RTVP1; CRISP7
- 储存: 冷冻(-20℃)
- 相关类别: 抗原
- 概 述:

Fusion protein corresponding to a region derived from 22-230 amino acids of human GLIPR1

技术规格:

Full name:	GLI pathogenesis-related 1
Synonyms:	GLIPR; RTVP1; CRISP7
Swissprot:	P48060
Gene Accession:	BC012510
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a protein with similarity to both the pathogenesis -related protein (PR) superfamily and the cysteine-rich secretory prot ein (CRISP) family. Increased expression of this gene is associated wi th myelomocytic differentiation in macrophage and decreased expres sion of this gene through gene methylation is associated with prosta te cancer. The protein has proapoptotic activities in prostate and bla dder cancer cells. This gene is a member of a cluster on chromosom



e 12 containing two other similar genes. Alternatively spliced variants which encode different protein isoforms have been described; howev er, not all variants have been fully characterized.